




**Welcome to
the Official Release of the
NIST23 Mass Spectral Libraries**

June 5, 2023

**Mass Spectrometry Data Center
National Institute of Standards and Technology
Gaithersburg, Maryland**


About the National Institute of Standards and Technology

- **1901: Founded as the National Bureau of Standards**
- **The national metrological institute (NMI) of the USA**
- **A non-regulatory agency of the US Department of Commerce**
- **Mission: To promote U.S. innovation and industrial competitiveness by advancing measurement science, standards, and technology in ways that enhance economic security and improve our quality of life**
- **Providing Standard Reference Data (SRD) one of the core functions of the Institute**
- **1968: Congress passed the Standard Reference Data Act**
- **2023: NIST currently provides more than 100 SRD products**



About the NIST Mass Spectrometry Data Center

- 1945: NIST produced the first mass spectral reference spectra
- 1988: the NIH/EPA mass spectral library was handed over to NIST
- Since that time NIST has grown the library over seven-fold
- 2023: 12th release of the NIST/NIH/EPA library
- All measurements, evaluation, software development done in house
- Total staff of 37
- Divided evenly between chemists and programmers



1988


The NIST/EPA/MSDC Mass Spectral Database, Personal Computer Versions 1.0 and 2.0

Distributor and price: U.S. National Institute of Standards and Technology, Office of Standard Reference Data, Gaithersburg, MD 20899, U.S.A. (Version 1 US\$ 750, Version 2 US\$ 975, upgrade from Version 1 to 2 US\$ 225).

Technical specifications:
Computer: Version 1: IBM XT, AT, PS/2 or compatibles with Hercules monochrome or CGA, EGA color monitor for optional graphics display. Version 2: IBM AT, PS/2 or compatibles with VGA, EGA, CGA color display for optional graphics.

Operating system: MS DOS
Minimum memory: 512 K (Version 1, 2), 640 K (Version 2 with all options).

Peripherals required: Hard disk (Version 1: 8 to 14 Mbytes, Version 2: 9 to 22 Mbytes). Regular Epson compatible printer for text only. For graphics printing HP Laserjet + or compatible printer.






Triennial Lift Off


Presenters

| | |
|---------------------------|----------------------------|
| Bill Wallace | Introduction |
| Steve Stein | Overview/Software |
| Tytus Mak | Selecting Compounds |
| Weihua Ji | EI Library |
| Yufang Zheng | Environmental/PFAS |
| Xiaoyu (Sara) Yang | Tandem Library |
| David Sparkman | Users |


NIST 23 Mass Spectral Libraries




NIST/EPA/NIH EI MS Library Evolution




GC/MS References



El Spectra over Time



Compounds in Tandem Library



How do we Select Compounds to Add to Our Libraries?

Tytus Mak

~~**What** new compounds are in NIST23?~~

How were new compounds selected for NIST23?

Why were these new compounds added to NIST23?



ZINC



FOODB



Substances
Added To Food



BIOCYC
Genome Database Collection

DRUGBANK



Hazardous
Substances
Databank



BitterDB
Institute of Biochemistry, Food Science
The Hebrew University of Jerusalem



hmdb
The Human Metabolome Database




EPA
PARIS III



NutriChem



ChEBI



npatlas



Blood
Exposome
Database



WIKIPEDIA
The Free Encyclopedia




PhytoLab

Automated Compound Selection Pipeline



Merge via
InChIKey

Meta-
Compendium
500K Structures



Coverage

Depends on Sample & Application

EI


| Citation | New | NIST23 |
|------------------------|------|--------|
| Wikipedia | 1570 | 6227 |
| EPA Tox | 2969 | 7117 |
| Food DB | 582 | 4273 |
| EU Contaminants | 6263 | 15149 |
| Protein Data Bank | 1194 | 4716 |
| Human Metabolome DB | 1992 | 9393 |
| PFAS | 161 | 749 |
| Adams (Essential Oils) | 2136 | 2136 |

Tandem

| Citation | New | NIST23 |
|---------------------|------|--------|
| Wikipedia | 1618 | 6424 |
| EPA Tox | 3181 | 8146 |
| Food DB | 602 | 4491 |
| EU Contaminants | 6553 | 15818 |
| Protein Data Bank | 1246 | 4945 |
| Human Metabolome DB | 2071 | 9686 |
| PFAS | 90 | 116 |
| Drug Bank | 1130 | 3991 |


Software

MS Search




Callable by External Programs

MS Interpreter EI or Tandem



AMDIS (gc/ms)



API & MSPepsearch

New Feature

Separate Names and Spectra - EI

Compound Name

Spectra
Replicates
Derivatives
Stereiosomers

The screenshot shows the NIST MS Search 3.0 interface. The search term 'BENZOICACID' is entered in the search bar. The results list includes various derivatives of benzoic acid. The 'spectrum' section is expanded, showing 'replicate 1' through 'replicate 4', and 'precursor + TMS' is selected. The mass spectrum plot shows relative intensity versus m/z, with major peaks at 105, 135, and 179. The chemical structure of Benzoic Acid, TMS derivative is shown as an inset. The 'Names' tab is selected at the bottom.

Search Results:

- Benzoic acid, 3,4-dimethoxy, n-butyl ester
- Benzoic 3-chloro-N-ethoxy-2,6-dimethoxybenzimidic anhydric
- Benzoic acid**
- Benzoic acid, 10,13-dimethyl-17-(1-methylpent-4-ynyl)-2,3,4,
- Benzoic acid, 10-chlorodecyl ester
- Benzoic acid, (1-(1,1'-biphenyl)-4-ylethylidene)hydrazide
- Benzoic acid, (1-[1,2-di(methoxymethoxy)-3-oxopropyl]-2-prc
- Benzoic acid, 1,1-dimethylethyl ester
- Benzoic acid, 1,2,3,4,5-pentamethylcyclopenta-2,4-dienyl est
- Benzoic acid 1,2,3,9-tetrahydropyrrolo[2,1-b]quinazolin-3-yl e
- Benzoic acid, [1-(2-methoxyethyl)-4-piperidyl] ester
- Benzoic acid, 1-(2-t-butyl-1-methyl-5-oxoimidazolidin-4-yl)-2-
- Benzoic acid, 1-(2-thienyl)ethylidenamino ester
- Benzoic acid, 1-(3,3-dimethyl-2-oxobicyclo[4.1.0]hept-7-yl)-e
- Benzoic acid, 1,3,3-trimethyl-2-oxa-bicyclo[2.2.2]oct-6-yl este
- Benzoic acid 1,3,3-trimethyl-4-phenyl-piperidin-4-yl ester
- Benzoic acid, 1,3-dioxan-5-yl ester


Mass Spectrum Data:

| m/z | Intensity |
|-----|-----------|
| 15 | 6 |
| 26 | 2 |
| 27 | 9 |
| 28 | 11 |
| 29 | 9 |
| 30 | 1 |
| 31 | 4 |
| 37 | 1 |
| 38 | 3 |
| 39 | 8 |
| 42 | 4 |
| 43 | 40 |
| 44 | 12 |
| 45 | 106 |
| 46 | 8 |
| 47 | 27 |
| 48 | 2 |
| 49 | 3 |
| 50 | 47 |
| 51 | 148 |
| 52 | 8 |
| 53 | 9 |
| 54 | 2 |
| 55 | 7 |
| 56 | 2 |
| 57 | 4 |
| 58 | 8 |
| 59 | 28 |
| 60 | 7 |
| 61 | 8 |
| 62 | 3 |
| 63 | 5 |
| 64 | 2 |
| 65 | 7 |
| 66 | 2 |
| 67 | 3 |
| 68 | 1 |
| 69 | 2 |
| 70 | 2 |
| 71 | 2 |
| 72 | 7 |
| 73 | 55 |
| 74 | 35 |
| 75 | 50 |
| 76 | 43 |
| 77 | 680 |
| 78 | 47 |
| 79 | 4 |
| 81 | 2 |
| 83 | 2 |
| 85 | 3 |
| 89 | 31 |
| 91 | 12 |
| 92 | 3 |
| 93 | 4 |

Separate Names and Spectra - Tandem

Compound
Name

Spectra
Energies
Fragmentation
Resolution
Charge
In-Source




New non-Spectral Data

(distinguishes *isospectral* compounds)

- **Reliable RI For All EI Library Compounds (40% => 100%)**
 - AIRI Fills Gaps (Lewis Geer)
 - RI Can Be Used for Scoring
- **'Other DB' Occurrences ('Prior Probability')**
 - Non-Spectral Databases Citing Each Compound
 - 61 DBs, Wide Applications (Environmental, Drug, Food, Metabolite, ...)
 - EI 21%, Tandem 78%

Retention Index Estimation - Old vs. New



New Search Type

Isomer ID by RI, Confirmation by 'Other DBs'


Search With o-Xylene Spectrum

No RI

| # | Lib | MF | RMF | Prob | RI | DBs | Name |
|-----|-----|-----|-----|------|--------|--------|-------------------------------|
| ⊕ 1 | M | 941 | 942 | 29.7 | 865 | 32 ... | p-Xylene |
| ⊕ 2 | M | 937 | 939 | 26.5 | 866 | 32 ... | Benzene, 1,3-dimethyl- |
| ⊕ 3 | M | 924 | 925 | 18.0 | 888 | 36 ... | <u>o-Xylene</u> |
| ⊕ 4 | M | 888 | 895 | 6.09 | 855 | 35 ... | Ethylbenzene |
| ⊕ 5 | M | 878 | 881 | 4.55 | 121... | | Benzeneethanol, α,β-dimethyl- |

Use RI

| # | Lib | MF | RMF | Prob | RI | origMF | DBs | Name |
|-----|-----|-----|-----|------|------|--------|--------|--------------------------------------|
| ⊕ 1 | M | 924 | 925 | 90.6 | 888 | 924 | 36 ... | <u>o-Xylene</u> |
| ⊕ 2 | M | 826 | 877 | 6.70 | 878* | 876 | | Cyclopentene, 1-ethenyl-3-methylene- |
| ⊕ 3 | M | 767 | 939 | 1.22 | 866 | 937 | 32 ... | Benzene, 1,3-dimethyl- |
| ⊕ 4 | M | 761 | 942 | 1.03 | 865 | 941 | 32 ... | p-Xylene |
| ⊕ 5 | M | 699 | 750 | 0.17 | 898 | 749 | 1 G | 1,3,5-Cyclooctatriene |



Hybrid Search

For compounds
not in the Library

The Hybrid Search: A Mass Spectral Library Search Method for Discovery of Modifications in Proteomics

Meghan C. Burke,^{*,†} Yuri A. Mirokhin,[†] Dmitrii V. Tchekhovskoi,[†] Sanford P. Markey,[†] Jenny Heidbrink Thompson,[‡] Christopher Larkin,[‡] and Stephen E. Stein[†]

Combining Fragment-Ion and Neutral-Loss Matching during Mass Spectral Library Searching: A New General Purpose Algorithm Applicable to Illicit Drug Identification

Arun S. Moorthy,^{*,†} William E. Wallace,[†] Anthony J. Kearsley,[‡] Dmitrii V. Tchekhovskoi,[†] and Stephen E. Stein[†]

Hybrid Search: A Method for Identifying Metabolites Absent from Tandem Mass Spectrometry Libraries

Brian T. Cooper,^{*,†,‡} Xinjian Yan,[‡] Yamil Simón-Manso,[‡] Dmitrii V. Tchekhovskoi,[‡] Yuri A. Mirokhin,[‡] and Stephen E. Stein[‡]

Increasing the Coverage of a Mass Spectral Library of Milk Oligosaccharides Using a Hybrid-Search-Based Bootstrapping Method and Milks from a Wide Variety of Mammals

Concepcion Africano Remoroza,^{*} Yuxue Liang, Tytus D. Mak, Yuri Mirokhin, Sergey L. Sheetlin, Xiaoyu Yang, Joice V. San Andres, Michael L. Power, and Stephen E. Stein

Structure Annotation of All Mass Spectra in Untargeted Metabolomics

Ivana Blaženović,[†] Tobias Kind,[†] Michael R. Sa,[†] Jian Ji,[‡] Arpana Vaniya,[†] Benjamin Wancewicz,[†] Bryan S. Roberts,[†] Hrvoje Torbašinović,[§] Tack Lee,^{||} Sajjan S. Mehta,[†] Megan R. Showalter,[†] Hosook Song,^{||} Jessica Kwok,[†] Dieter Jahn,^{⊥,#} Jayoung Kim,^{∇,○,◆,¶} and Oliver Fiehn^{*,†}

LC-MS/MS Software for Screening Unknown Erectile Dysfunction Drugs and Analogues: Artificial Neural Network Classification, Peak-Count Scoring, Simple Similarity Search, and Hybrid Similarity Search Algorithms


Inae Jang,[†] Jae-ung Lee,[†] Jung-min Lee,[†] Beom Hee Kim,[‡] Bongjin Moon,[†] Jongki Hong,^{*,‡}

NIST

Non-NIST

Hybrid Search

28 Da (CH₂CH₂) Difference



Formula Calculator

Mass + = 182.1177 40 ppm Loss = 122.0366


Formula C17H21NO4(H) Parent = 304.1543

| 1 Ion | O+E | RDB | Mass | Ppm | C | H | N | O |
|-----------|------|-----|-----------|------|----|----|---|---|
| C10H16NO2 | Even | 3.5 | 182.11756 | -0.8 | 10 | 16 | 1 | 2 |

Formula/Mass/Isotope Calculator

Maximum Rate: 101(chain) @ 182.1176 m/z

| m/z | exact m/z | formula | loss | type | rate | rel.rate | ion.rate | ppm | abund | O+E |
|-----------|------------|-----------|--------|----------------|------|----------|----------|------|-------|------|
| 182 (1/2) | 182.117556 | C10H16NO2 | C7H6O2 | H-Displacement | 101 | 88 | 79 | -0.8 | 999 | Even |




Fragmentation Class Rates/Formulas/...

Structure for Selected Peak Appears in Red

MS Interpreter
High or Low Res
Called from NISTMS


Mass Spectrum for C17H21NO4(H)⁺ Mass = 304.1543; CAS = 50-36-2; [M+H]⁺ protonated Cocaine [M+H]⁺





NIST/EPA/NIH EI Library: Evaluation

Weihua Ji


Expanding the EI MS Library





Evaluating Spectra with MS Interpreter





RI and AIRI Assist Spectrum Evaluation



Experimental RI: 3027




Evaluating Spectra with the Hybrid Search



PFAS & Environmental


Yufang Zheng

PFAS Selection




NIST MS Interpreter Aids Spectrum Evaluation

This program determines the probable origin of peaks using thermochemical rules




Name: 4,4,5,5,5-Pentafluoro-1-pentanol, TMS
Formula: C₈H₁₅F₅OSi
MW: 250




Hybrid Search to Aid Spectrum Evaluation


4,4,5,5,5-Pentafluoro-Pentanol, TMS



Hybrid search spectrum



4,4,5,5,6,6,6-Heptafluoro-hexanol, TMS





**400 PFAS spectra for 231 PFAS
Added to NIST23**


6/7/23 WOE 2:50PM GC/MS: Instrumentation and Applications; Room 332

*Enhancing the Coverage and Quality of Spectra of per- and polyfluoroalkyl
Substances in a Comprehensive Electron Ionization Mass Spectral Library*



NIST23 Tandem Mass Spectral Library

Xiaoyu (Sara) Yang

NIST23 Tandem Mass Spectral Library



- Features: **high quality**
comprehensive
high and low resolution
library searching software
- Application: metabolomics, health, environmental, pharmaceutical, forensics, food, and agriculture, etc.



Identifying Compounds by Searching the *NIST23* Tandem MS Library

NIST MS Search 3.0 (Tandem mode) - [MS/MS, Presearch Default - 97 spectra]

File Search View Tools Options Window Help

MS Search 3.0

1. Cluster_017677

| # | Src. | Name |
|---|------|----------------|
| 1 | A | Cluster_017677 |
| 2 | A | Cluster_017678 |
| 3 | A | Cluster_017679 |

hr_msms_nist; hr_msms_nist#2; 1934658 total spectra

(Text File) Cluster_017677

Name: Cluster_017677
MW: N/A ID#: 13 DB: Text File
Comment: Spectrum=Consensus Parent=253.2161 Nreps=6/6 RT=22.4 Ion_mode: P
Collision_energy: 20% ID: 17677 Sample: plasma_hcd_pos
209 m/z Values and Intensities:
55.0178 1.94 | 55.0544 174.11 | 56.0496 1.21 | 57.0335 21.85 |

Human Plasma

NIST23 Library

ω-6-Hexadecenlactone

| # | Lib. | Score | DotProd | Rev-Dot | Prob. |
|---|------|-------|---------|---------|-------|
| 1 | hr | 904 | 916 | 935 | 73.8 |
| 2 | hr | 867 | 880 | 933 | 24.4 |
| 3 | hr | 724 | 772 | 831 | 0.80 |
| 4 | hr | 703 | 752 | 924 | 0.42 |
| 5 | hr | 696 | 733 | 901 | 0.34 |
| 6 | hr | 660 | 734 | 873 | 0.11 |
| 7 | hr | 635 | 770 | 876 | 0.05 |
| 8 | hr | 616 | 663 | 767 | 0.03 |
| 9 | hr | 612 | 676 | 808 | 0.02 |


(hr_msms_nist#2) ω-6-Hexadecenlactone [M+H]⁺ HCD 20% P=253.2

Name: ω-6-Hexadecenlactone
Precursor type: [M+H]⁺
Instrument type: HCD
Collision energy: NCE=20% 19eV
Precursor m/z: 253.2162
Formula: C₁₆H₂₈O₂
MW: 252 Exact Mass: 252.20893 CAS#: 7779-50-2 NIST#: 4419512 ID#: 611023 DB: hr_msms_nist#2
Comment: NIST Mass Spectrometry Data Center
Related CAS#: 223104107
Ion mode: P
Instrument: Orbitrap Fusion Lumos
Ionization: ESI
Collision gas: N₂
Sample inlet: direct flow injection
Spectrum type: MS2





Lib. Search Other Search Names Compare Librarian

Type of Search: MS/MS Displayed: MS/MS

How Did We Create Our High Quality and Comprehensive Library?



561 Extractable & Leachable Compounds with APCI and 352 Compounds with APCI and ESI



NIST End-User and Distributor Support

Contractor

O. David Sparkman

massspec@nist.gov

Or



ods@csi.com

1-209-483-5740

Direct Benefit to Distributors

- **Compare NIST 23 with older versions and benefits of upgrading**
- **Interface MS Search to 3rd party software**
- **Use MS Search to perform tasks done by other software**
- **Clarify NIST MS Data licensing requirements**
- **Technical support for MS Search for your end users**
- **Submit inquires to massspec@nist.gov**

Custom Branding For Distributors




Branding Accomplished by:

Adding text file with name **Serial.dll**
to installation media

Contact me for details

NIST 23 & MS Search Manual



Ver30Man.PDF
In
C:\NIST23\MSSearch\

Windows Start Menu

NIST Standard Reference Database 1A

Mass Spectral Libraries (NIST/EPA/NIH EI & NIST Tandem 2023) and NIST Mass Spectral Search Program (Version 3.0)

For Use with Microsoft® Windows

User's Guide

The NIST Mass Spectrometry Data Center

Stephen E. Stein, Technical Leader

William E. Wallace, Group Leader

Evaluation:

Weihua Ji (EI)

James L. Little (EI)

H. Martin Garraffo (Tandem)

Sanford P. Markey (EI)

W. Gary Mallard (EI)

Pedatsur Neta (Tandem)

Kirill V. Tretyakov (EI)

O. David Sparkman (EI)

Measurement/Evaluation:

N. Rabe Andriamaharavo (EI)

Tallat H. Bukhari (Tandem)

Edward P. Erisman (EI)

Yuxue Liang (Tandem)

Yi Liu (Tandem)

Yamil Simòn-Manso (Tandem)

Quan-Long Pu (EI)

Concepcion A. Remoroza (Tandem)

Nino G. Todua (EI)

Yufang Zhang (EI)

Hyperlinks in Manual and to External files

Hyperlinks

And
Delta_Mass_Guide.PDF
Delta Mass Guide for EI

The Group/Element agreements below have been observed in the evolution of actual spectra using the Hybrid Search

There are **Group/Element (1)** and **Group Element (2)**. If the sign of δ is + (positive, assumed) then **Group Element (2)** is a possible replacement for what may be **Group/Element (1)** in the structure of the library compound used as the model. If the sign of δ is - (negative, specified) then Group Element (1) is a possible replacement for what may be Group/Element (2) in the structure of the library compound used as the model

Hyperlink in Mar OtherDBs.PDF

| Catalog Name |
|--------------------------------------------------|
| EPACOMPTOX_EPA official InertIngredients_Nov2019 |

| |
|---------------------------------|
| EPACOMPTOX_Pesticides |
| InertFinder_Nov2020 |
| EPACOMPTOX_PLASTICS NORMAN 2019 |
| Thermo_AdditivesList2022 |

DeltaMass Nominal

0
1
1
1
1
1
1
1
1
1
1
1

imidazole ring
 CH3SO2 group
 methyl on aromatic ring
 CH2NH2 group
 nitrogen in heterocyclic aromatic ring
 amine on aromatic ring
 N in six membered heterocyclic aromatic
 insertion of N in place of carbon in 5-membered heterocyclic ring
 TBDMS derivative attached to two NH groups
 CH2NH2 onaromatic ring

Group/Element (1)

Group/Element (2)

furan ring
 NH2SO2 group
 amine on aromatic ring
 CH2CH3 group
 phenyl aromatic ring
 phenol ring
 aromatic ring no nitrogen incorporated
 no insertion
 TDDMS derivative inserted on O and one in group on aromatic ring
 CH2OH on aromatic ring

Contents

| | |
|----------------------------------------|----|
| Preface | 1 |
| Nominal Mass for Libraries | 3 |
| OS Compatibilities | 5 |
| Installation | 6 |
| FIRST QUESTION | 8 |
| Getting Started | 9 |
| First Action After Program Setup | 10 |
| | 11 |
| | 12 |

URL

c.epa.gov/dashboard/chemical-lists/PESTINERTS

Both in:

c.epa.gov/dashboard/chemical-lists/PESTINERTS
C:\NIST23\MSSearch


c.epa.gov/dashboard/chemical-lists/PPDBLISTB

Training Videos for NIST 23 and beyond

<https://littlesandsailing.wordpress.com/>

A "Little" Mass Spec and Sailing
Organic Mass Spectrometry, NMR, Sailing, Tesla, Duplicate Bridge

Mass Spec Interpretation Services



The graphic features a mass spectrum plot on the left with peaks at m/z 15, 16, 27, 31, 39, 42, 44, 57, and 58. A 'MASSMAN' license plate is shown above the plot. To the right is a diagram of a mass spectrometer with a magnifying glass over the m/z 569.2876 peak. A computer monitor displays 'CAS NO. 20583-873'. A chemical structure of acetone is also present. The text 'C₆H₁₂O₆ and CASE COATING' is written near the magnifying glass.

About Me | My Topics | Others Links

Overview "Known Unknowns" | Posted by: [tvasailor](#) | May 24, 2012

4/27/23: NIST2023 Libraries and Software Available in Early June
4/8/2023: FREE Drawing Program, "Mass Spec Calculator Professional"

1/25/2023: "Approaches for Identifying New Psychoactive Substances (NPS)" at 2023 Current Trends in Seized Drugs Analysis Symposium

11/24/22: Creating and Using Retention Indices in NIST Software
11/15/22: Forensic@NIST2022 Workshop: Seized Drug Analyses
11/01/22: Training for Wiley KnowItAll New EI Mass Spec Capabilities
10/01/22: FREE Courses for Unknown Identification Using NIST Search

BROWSE

Monthly Archives

Search

LINKS

- Accurate Mass NIST User Libraries
- Archives of Early MS
- Bridge Card Game

Training videos, handouts and instruction on NIST Libraries and NIST Mass Spectral Search Program

See You at Booth 713



- Speaker Order
 - Bill (2), Steve (7), Tytus (3), Steve (10-1,1,2,3,2,1)
 - Weihua (4), Yufang (4), Sara (4), Sparkman (5)
- Stein: Coverage, Software, Names/Spectra, Search Methods, Non-Spectral Data, EI-RI, Search Types, Hybrid Search, MS Interpreter