


Single Step Protein Extraction from Trace Amount Human Hair for Genetically Variant Peptide Detection

Zheng Zhang, Meghan C. Burke, William E. Wallace, Yuxue Liang, Sergey L. Sheetlin, Yuri A. Mirokhin, Dmitrii V. Tchekhovskoi, Stephen E. Stein
National Institute of Standards and Technology, Mass Spectrometry Data Center, Gaithersburg, Maryland, United States

Introduction


Recent reports have demonstrated that genetically variant peptides (GVPs) derived from human hair proteins can be used to differentiate individuals of different biogeographic origin (PMID:27603779 DOI:10.1371/journal.pone.0160653). We report a single-step direct extraction method ('Direct' Method) for hair proteins that is significantly more sensitive than two previously published methods (PMID:27603779 DOI:10.1371/journal.pone.0160653 and PMID:27741315 DOI:10.1371/journal.pone.0164993). It is also simple to implement and gives reproducible results. Together with the construction of a hair specific peptide mass-spectral library including previously reported GVPs, it provides a means of reliable analyses of human hair proteome which can utilize trace amount forensic samples.

Result – Test A. Comparison of Reproducibility




Test A: The two gel images compare the reproducibility of single-step Direct method and NaOH+SDS method using 5 cm-long hair shaft samples from the same individual donor across 8 replicates (Left: A to H; Right: 1A to 1H). A MW standard was loaded in the first lane. Note that the NaOH+SDS gel includes a 9th lane for which the extraction from ten 5cm-long hair shaft samples was included as a reference. The major bands that correspond to type I and type II hair cuticular keratins were labeled.

Result – Test B. Comparison of Artefacts



Test B: Comparison of experimentally introduced artifactual modifications among three methods using our recently developed hybrid search (PMID:28367633 DOI:10.1021/acs.jproteome.6b00988): Cleavable Surfactant method (red), NaOH+SDS method (green) and single-step Direct method (blue). The compared experimentally introduced artifactual modifications chosen as examples are: acetaldehyde (upper left), acetylation (upper right), formylation (lower left) and over alkylation (lower right).

Result – Test D. Reduction of Hair Length




Test D: This gel image (left) shows the separation of hair proteins from 5, 2.5, and 1 cm-long hair shaft samples from the same individual donor. A MW standard was loaded in the first lane. Bands for type I and type II hair cuticular keratins were labeled. The table (right) summarizes protein, peptide, and GVP ion identifications from 5, 2.5, and 1 cm-long hair shaft samples.

Methods


- Samples were commercially obtained. Five cm-long hair shaft sample was combined with 50 μ l of the commercial sample loading buffer and 50 mM dithiothreitol (DTT). After heating at 90°C for 30 min, extracted hair proteins were loaded onto a protein gel and then separated by size. We call this single-step protein extraction by heating a single hair shaft in sample loading buffer the 'Direct' method. After gel staining, in-gel digestion was performed and extracted peptides were injected to a mass spectrometer for LC-MS/MS analyses.

Single-Step Direct Extraction Method



- We compared this single-step Direct extraction method to two previously published methods: 1) NaOH with SDS repeated extraction ('NaOH + SDS') method (PMID:27741315 DOI:10.1371/journal.pone.0164993); 2) Cleavable Surfactant (ProteaseMax) based method (PMID:27603779 DOI:10.1371/journal.pone.0160653)
- Using the raw data files generated from the study, a hair specific peptide mass-spectral library was constructed. It has a sequence coverage of hair cuticular keratins of about 70% and includes 14 previously reported GVPs.

Result – Test C. Comparison of Sensitivity



Test C: The sensitivity of all three methods were measured by comparing multiple metrics across a dilution series from 5D to 1280D: (upper left) the total number of ions; (upper right) the total number of peptides; (lower left) total number of proteins; (lower right) total number of published GVP ions detected in mass spectral data from 5 cm-long hair shaft sample derived proteins that were extracted using the Cleavable Surfactant method (red), NaOH+SDS method (green) or the Direct method (blue). Actual data has been labeled on the points of each dilution series.

Summary & Conclusions

- Our results show that single-step Direct method is the most sensitive and reliable method among the three tested methods.
- This finding was confirmed by determining that differences in the results of three tested methods are not due to experimentally introduced artefacts including modifications.
- The construction of a human hair specific peptide mass spectral library including previously reported GVPs provides a means of reliable identification of the human hair proteome, adding another layer of sensitivity to our method, for quick and accurate peptide identification and GVP detection.