


chemdata.nist.gov


Introduction

- Previously¹, multiple tools developed at the NIST Mass Spectrometry Data Center were integrated into a single graphical user interface (GUI)
- This allowed creating series of sequentially operating programs to be configured and executed through a pre-defined 'pipeline'
- This poster describes extension of this pipeline for the specific task of comprehensive glycopeptide analysis and building high-quality glycopeptide mass spectral libraries from raw data
- The pipeline is implemented as a GUI with minimum required input parameters
- The generated libraries can be used with an updated version (not currently available for download) of the NIST Search Software²
- The NIST Search Software can be used for examining and comparing spectra in the resulting glycopeptide libraries




Types of generated libraries

(1) MS/MS glycopeptide libraries; (2) Glycopeptide abundance distribution spectra (GADS³) libraries; (3) Glycan libraries; (4) MS/MS peptide libraries.



Creating a data analysis pipeline for producing high-quality glycopeptide mass spectral libraries

Sergey L. Sheetlin, Yuri A. Mirokhin, Guanghui Wang, Xiaoyu Yang, Concepcion A. Remoroza, Yi Liu, Xinjian Yan, Dmitrii V. Tchekhovskoi, Zachary C. Goecker, Meghan C. Burke, Stephen E. Stein Biomolecular Measurement Division, National Institute of Standards and Technology (NIST), Gaithersburg, MD 20899, USA



- Computation can be started from any step using the button "Run"
- The pipeline GUI analyzes the input settings and outputs warnings if some problems are detected during execution
- A batch software for processing multiple input data is available

identified are necessarily the best available for the purpose.

| projects | | |
|--|--|---------------------|
| | Additional parameters | omputation status |
| | Start calculation from this step | omputation status |
| folder Select | Prev. Cur. | omputation can be |
| databases | | arted from any step |
| g-483-1997).fasta | Formatting MSFragger output | |
| pen Remove | XICAnalyzer | |
| | Validating spectra | |
| xt | Making GADS | |
| | Building glycan library | Detailed status for |
| plete MSMS library | O Building NIST MS/MS libraries | each file and step |
| , | Details | |
| ut file | Compare parameters with defaults | |
| ascot converter | | |
| | Custom enzymes | Enzyme 1 Enzyme 2 |
| -10-25_380-2000_120_HS -09-17_380-2000_120_HS | Enzyme's cutting terminal (search_enzyme_sense) | C ~ C ~ |
| -09-23_380-2000_120_H5 -09-17_380-2000_120_H5 | Name of enzyme to be written to the pepXML file (search_enzyme_name) | trypsin_gluc |
| | Residues after which the enzyme cuts (search_enzyme_cut) | DEKR |
| ~ | Residues that the enzyme will not cut before (search_enzyme_nocut) | Р |
| Show relative paths | Allowed number of missed cleavages per peptide (allowed_missed_cleavage) | 3 🔺 |
| hoose protease | | |
| | | |
| | Screen output | of |
| | programs | |
| ut Run | Refresh controls | |
| | | .: |
| | | |
| Doguin | omonto | |

κequirements

References

(¹) Sheetlin, Sergey L., at al, Graphical User Interface for Integration of NIST Mass Spectrometry Data Center *Library Building Tools.* Poster on ASMS 2021 Conference on Mass Spectrometry and Allied Topics.

(³) Remoroza, C.A., et al., *Representing and Comparing Site-Specific Glycan Abundance Distributions of*

(⁴) Kong, A. T., et al., *MSFragger: ultrafast and comprehensive peptide identification in mass spectrometry*-

(⁶) Rudnick, P.A., et al., A Description of the Clinical Proteomic Tumor Analysis Consortium (CPTAC) Common

(⁷) Yang, X., et al., MS_Piano: A Software Tool for Annotating Peaks in CID Tandem Mass Spectra of Peptides and

(⁸) Meghan C. Burke, et al., The Hybrid Search: A Mass Spectral Library Search Method for Discovery of

Restrictions and Disclaimers

Certain commercial equipment, instruments, or materials are identified in this document. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the products